## **SAFETY RECALL**

## **Drive Recalls Bed Assist Handles and Rail**

**Hazard:** Users can become entrapped within the bed rail or between the rail and the side of the mattress. This poses a serious entrapment hazard and risk of death by asphyxia.



RTL15063-ADJ (Home Bed Assist Handle)



15064 (Bed Assist Handle)



RTL15062 (Bed Assist Rail w/Folding Board)



RTL15073 (Home Bed Assist Handle)

**Product Description:** The recall involves 3 models of bed assist handles and 1 bed assist rail which were sold between 2007 and December 2021.

**Remedy:** Consumers should immediately stop using the recalled products and contact Drive to receive a refund.

## **For More Information**

Visit <u>www.recallrtr.com/bedhandles</u> or call 1-877-467-3099 from 8 a.m. to 7 p.m. ET, Monday through Friday.